



USL16-NGL: Operational Throughput Improvements

**Michael Foxe, Theodore Bowyer,
Ian Cameron, Matthew Cooper,
James Hayes, Michael Mayer,
Jennifer Mendez, Johnathan Slack**

CTBTO Science and Technology Conference
Vienna, Austria, June 24-28, 2019



PNNL is operated by Battelle for the U.S. Department of Energy

Approved for public release; distribution is unlimited

The authors wish to acknowledge the funding support of the Nuclear Arms Control Technology Program of the U.S. Department of Defense, Defense Threat Reduction Agency.


The views expressed here do not necessarily reflect the opinion of the United States Government, the United States Department of Energy, or the Pacific Northwest National Laboratory

Outline

- USL16-NGL
 - Current System Overview and Process
- Sample Throughput Improvements
 - Motivation
 - Methods
- Detector Calibration and Analysis


USL16-NGL

- Security
 - Certified Personnel
 - Chain of Custody
- Environmental Monitoring
 - Temperature & Humidity
 - Radon – Not needed for beta-gamma (see it if it is present)
- Documentation
 - Quality Control
 - Serial Numbers and Calibration Certificates
 - Maintenance Log
 - Only make approved updates
- Software/Calculation
 - Version Control
 - Automated sample processing



Current USL16-NGL Detector

- Beta gamma system
 - Plastic beta cell
 - NaI
- Dual PMT readout for both the plastic scintillator and NaI
- Gain matching of PMTs performed during calibration process





Motivation for Higher Throughput

- Multiple samples in a Proficiency Test Exercise (PTE)
- Multiple station samples shipped at once
- Memory effect from calibration or PTE spikes
- Two solutions:
 - Coated cells
 - ✓ Solves the memory effect issues, but still requires gas backgrounds during sample processing
 - Additional detectors
 - ✓ See memory effect if not coated, but can alternate high and low-activity detectors
 - ✓ Can convert serial sample processing to parallel
 - ✓ Allows for additional re-verification of a measurement if there is uncertainty on the results



Implementation of Four Additional Detectors

- Single PMT readout instead of double ended beta cell
 - Fewer PMTs
 - Improved energy resolution
 - Smaller detector volume



“Quad” Detector Setup

- Small volume manifold and cryotrap
- Additional roughing pump to purify samples




Required System Improvements

- Slightly smaller cells
 - Must be accounted for when performing volumetric transfers
 - ✓ Required a new lower-manifold for the additional detectors
- Plumbing connection between the processing manifold and the lower quad detector manifold
 - Use nitrogen to push the xenon down to the lower manifold
- Retain the recollection ability if we want to re-count the sample
 - Secondary recollection port to minimize sample loss

Gas Manifold

- ▶ Xenon volume is quantified using MKS 910 pressure sensors
- ▶ Xenon is collected on a liquid nitrogen cooled low volume, high surface area trap
- ▶ Xenon is pushed with nitrogen to the lower manifold
- ▶ Xenon is volumetrically transferred to a beta-gamma detector



Sample Processing Upgrades

	00:00	02:00	04:00	06:00	08:00	10:00	12:00	14:00
Sample processing		Process sample #1	Process sample #2	Process sample #3	Process sample #4	Process sample #5		
ARSA cell	Gas Background	QC Count	Count Sample #1					
Quad cell 0	Gas Background		QC Count	Count sample #2				
Quad cell 1	Gas Background		QC Count	Gas Background	Count sample #3			
Quad cell 2	Gas Background		QC Count	Gas Background		Count sample #4		
Quad cell 3	Gas Background		QC Count	Gas Background			Count sample #5	


Control Page – Automation of Processes

Current State	Last Status Message	GUI Lock	Send Mail
purify-sample	Waiting for LNLevel to come up to 4.00 inches. Currently at: 0.03	04:17	

[Control Page](#)
[Flow Schematic](#)
[Status](#)
[Detector](#)
[Realtime Plotter](#)
[Event Viewer](#)
[SoH Viewer](#)
[System Volumes](#)
[About RL16](#)

Parameter	Value	Default	Units	Last Modified
COLLECTION_END_PRESSURE	2.18	2	torr	2016/06/20-15:19:18
COLLECTION_PREP_INTERVAL	30	30	s	
COLLECTION_START_INTERVAL	30	30	s	
COLLECT_RESIDUAL_SAMPLE_INTERVAL	180	180	torr	2016/04/01-13:18:04
COLLECT_RESIDUAL_TRAP_TEMP	-170	-170	C	2016/04/01-10:26:54
ELUTION_EQILIBRATION_TIME	180	180	s	
ELUTION_TRAP_TEMP	27.0	27.0	C	2016/04/01-13:18:04
FINAL_PNF_CYCLE_COUNT	3	3	none	2016/06/20-13:50:08
FINAL_PNF_FINAL_PUMP_INTERVAL	600	600	s	2016/06/16-09:37:17
FINAL_PNF_FLUSH_INTERVAL	60	60	s	2015/06/15-14:35:44
FINAL_PNF_PUMP_INTERVAL	600	600	s	2015/05/21-13:55:07
GAS_BACKGROUND_INTERVAL	79200	86400	s	2016/06/21-16:47:06
HISTOGRAM_COUNT_INTERVAL	86400	86400	s	2016/06/08-13:53:22
HISTOGRAM_PRELIM_INTERVAL	7200	7200	s	2016/04/01-13:18:04
INITIAL_PNF_CYCLE_COUNT	3	3	none	2016/06/03-11:28:20
INITIAL_PNF_FINAL_PUMP_INTERVAL	120	120	s	2016/06/03-11:28:20
INITIAL_PNF_FLUSH_INTERVAL	10	10	s	2015/05/21-13:55:07
INITIAL_PNF_PUMP_INTERVAL	120	120	s	2015/02/17-17:17:36
LOAD_CELL_FLUSH_INTERVAL	15	15	s	2015/06/17-15:25:49
LOAD_CELL_VOLUMETRIC_INTERVAL	120	120	s	2015/06/17-15:53:50
LOAD_CELL_VOL_FLUSH_PREP_PRESSURE	250	250	torr	
PUMP_AND_FLUSH_PRESSURE	300	300	torr	2016/02/25-14:49:56
PURIFY_CYCLE_COUNT	3	3	cycles	2016/06/20-13:50:08
PURIFY_FREEZE_INTERVAL	180	180	s	2016/04/01-13:18:04
PURIFY_PUMP_INTERVAL	45	45	s	2013/09/23-15:56:14
PURIFY_THAW_INTERVAL	180	180	s	2016/06/03-11:19:04
QC_COUNT_INTERVAL	1200	1200	s	2016/06/03-11:28:20
REGEN_CARRIER_INTERVAL	20	20	s	
REGEN_CARRIER_PRESSURE	100	100	torr	2016/02/29-14:27:01
REGEN_COOL_DOWN_TEMP	100	100	C	
REGEN_TRAP_TEMP	100	100	C	
REGEN_VACUUM_INTERVAL	600	600	s	2016/02/26-12:07:24
TRAP_COOL_DOWN_INTERVAL	300	300	s	2015/05/21-13:55:07
TRAP_LN_LEVEL	4.0	4.0	in	2016/06/08-13:53:22
UNLOAD_QUANTIFICATION_INTERVAL	30	30	s	2015/06/22-09:36:06
UNLOAD_SAMPLE_FREEZE_INTERVAL	120	120	s	2015/06/22-09:36:06
UNLOAD_SAMPLE_PUMP_INTERVAL	45	45	s	2015/06/22-09:36:06

[Cancel Edits](#)
[Apply Edits](#)
[Revert All To Defaults](#)



Detector Utilization

- 4 new detectors – alternate calibration cycles
- 2 detectors are “high activity/regular duration” counts
- 2 are kept for “low activity/long duration” counts
- Which 2 are which can be alternated once there has been a period of no activity so they are all “low activity” detectors
 - Avoid any biasing between detectors

USL16_005		USL16_005	PTE	USL16_005	Calibration	USL16_005	PTE	USL16_005	Station	USL16_005
USL16_006	Calibration	USL16_006	PTE	USL16_006		USL16_006	Station	USL16_006	Station	USL16_006
USL16_007	Calibration	USL16_007	PTE	USL16_007		USL16_007	Station	USL16_007	Station	USL16_007
USL16_008		USL16_008	Station	USL16_008	Calibration	USL16_008	PTE	USL16_008	Station	USL16_008
USL16_009		USL16_009	Station	USL16_009	Calibration	USL16_009	PTE	USL16_009	Station	USL16_009

Memory Effect

- Consistent memory effect measurements between detector types

Table 6. Memory Effect Results: 24 hour measurements calculated for the same reference time








**Original Double
Ended Detector**

	Sample Activity Concentration (mBq/m ³)	Gas Background Activity Concentration (mBq/m ³)	Memory Effect (%)
Xe-133 (80)	931	25.9	2.78
Xe-133 (30)	909	22.6	2.48
Xe-133 (sum)	924	25.1	2.72







**Additional Single
Ended Detector**

Xe-131m	14.5	0.32	2.21
---------	------	------	------

Calibration – 5 Isotopes



Energy Resolution Improvement



Conclusions

- As the noble gas laboratories move to a higher utilization rate, increased sample throughput is desired
- USL16-NGL is implementing 4 additional detectors to improve the sample throughput
- The additional detectors have improved energy resolution, but no memory effect coating
- USL16-NGL will have the capability of processing 2-4 samples per day with the additional detectors

Thank you